



Email : sales@laboratoryinstrumentindia.com

Product Name :
PLC Programmable Logic Controller And Automation
Training System

Product Code :
LBNY-0005-12600015



www.laboratoryinstrumentindia.com

Description :

The PLC Programmable Logic Controller And Automation Training System has a built-in powerful industrial PLC employing the latest technology features which allows the management of all the typical variables of an industrial process. The PLC Programmable Logic Controller and Automation Training System is today the core of most of the control systems. The PLC programs can be managed by a PC and dedicated software for Windows. The system is supplied with manuals which describe all parts of the unit, the installation and utilization procedures, as well as many exercises with the relative results. In fact, very often in any automated system the commands and controls are processed and executed by PLC. The unit makes available in a single structure all the equipment required to carry out exercises at any level of complexity, from the simplest to the most demanding.

Technical Specification :

PLC Programmable Logic Controller And Automation Training System



Laboratory Instrument India