




Email : sales@laboratoryinstrumentindia.com

Product Name :
Communication Trainer

Product Code :
LBNY-0002-610004



Description :

Communication Trainer

Technical Specification :

Data Communication Trainer, the switch faults are provided on-board to study different effects on circuit and 20 Nos. of test points are provided on board to observe various intermediate signals.

Serial Communication : RS-232 port (9 pin & 25 pin)

Parallel Communication : 25 Pin LPT Port Wireless Communication : Infrared Transmitter : IR LED 920mm

Infrared Receiver : Direct TTL output

Baud Rate : 2.4 kbps (max)

Carrier Frequency : 38Khz

Communication Trainer, Communication Trainer Equipments, Communication Trainer Tools, Communication Trainer Tool Kits, Communication Trainer Manufacturers, Communication Trainer Suppliers from India, China, Kenya



Laboratory Instrument India

