



Email : sales@laboratoryinstrumentindia.com

Product Name :
Adequate Protractor For Fixing To The Track

Product Code :
LBNY-0004-780011



Description :

Adequate Protractor For Fixing To The Track

Technical Specification :

Mechanical bench accessory 09:

Adequate protractor for fixing to the track

The protractor has two sets of measurements. The marking on the inner circle is 0° to 180° in the opposite direction to the markings on the outer circle.

It is an instrument for measuring or constructing an angle of a given measure. It is a circular or semicircular piece of metal or plastic.



Laboratory Instrument India