



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Adequate Protractor For Fixing To The Track

**Product Code :**  
LBNY-0004-780011



**Description :**

Adequate Protractor For Fixing To The Track

**Technical Specification :**

Mechanical bench accessory 09:

Adequate protractor for fixing to the track

The protractor has two sets of measurements. The marking on the inner circle is 0° to 180° in the opposite direction to the markings on the outer circle.

It is an instrument for measuring or constructing an angle of a given measure. It is a circular or semicircular piece of metal or plastic.



**Laboratory Instrument India**