



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pressure Transducers Study Module

**Product Code :**  
LBNY-0005-12600017



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The automatic data acquisition and analysis system for Personal Computer allows to take the measurements automatically, plot on screen or print out and file the test data. The Pressure Transducers Study Module allows to study the calibration methods for a signal conditioner instrument connected to a pressure transducer and to verify the system linearity, repeatability and hysteresis.

**Technical Specification :**

Absorbed power: 0.3 kW max  
Electric supply: 230 V, 50/60 Hz, single phase.  
Weight and dimensions:  
Dimensions: 520 x 310 x 540 h mm  
Net weight: 13 kg.



**Laboratory Instrument India**