



Email : sales@laboratoryinstrumentindia.com

Product Name :
Acceleration Module

Product Code :
LBNY-0005-12600018



www.laboratoryinstrumentindia.com

Description :

The Acceleration Module is suitable to study the behaviour of a system by using a piezoresistive acceleration transducer.

Technical Specification :

Required services:

Electrical supply: 220/240V, 50/60 Hz

Weight and Dimensions:

Dimensions: 600x400x400 h mm

Weight: 30 kg

Experiments:

Acceleration Module is possible to study the system both under free and forced oscillation conditions in order to determine natural frequency and damping coefficient and compatibles or by oscilloscope.



Laboratory Instrument India

