



Email : sales@laboratoryinstrumentindia.com

Product Name :
Temperature Transducers Study Module

Product Code :
LBNY-0005-12600025



www.laboratoryinstrumentindia.com

Description :

A special software program operating in MS-Windows environment enable to perform the execution of tests through a Personal Computer. The Temperature Transducers Study Module allows the calibration of instruments connected to temperature transducers, such as thermo couples and thermo resistances, and it demonstrates experimentally the concept of measure uncertainty. The unit includes all the measuring instruments necessary for the execution of the tests and numerous measuring points to check the most important signals with the aid of external instruments.

Technical Specification :

Required services:

Electrical supply: 220 V, 50/60 Hz, single-phase, 2 kW

Dimensions: 480 x 410 x 600 h mm

Net weight: 32 kg.



Laboratory Instrument India