



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Modular Process Control Trainer

**Product Code :**  
LBNY-0005-12600026



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The basic module includes the electrical equipment and the components for the Level control. The Modular Process Control Trainer system is a modular and low cost trainer based on Personal Computer and suitable to study the Process Control techniques. It is possible to carry out the laboratory activity in an extremely productive way using the potential of a Personal Computer. Some additional kits allow to fulfill the Flow, Pressure, Temperature, pH and conductivity control by using most of the basic module components.

**Technical Specification :**

**Additional kits**

They are made by different elements with contained weight and dimensions.

Optional electronic regulator

Dimensions: 200x350x200 mm

Weight: 1.5 kg

**Weight and Dimensions:**

Base system

Dimensions: 700x400x1000 mm

Weight: 50 kg approx.

---

# Laboratory Instrument India

