



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Stepping Motor Control Module

**Product Code :**  
LBNY-0005-12600027



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Stepping motor, 200 step/rev, max torque = 60 N· cm  
No. 2 limit micro switches  
Graduated disc, 0 360  
Encoder (200 pulse/rev)  
Gear motor 90/1 rev.  
Electric supply box consisting of:  
a) Ammeter for measuring the motor absorbing current  
b) Diagram on front panel showing the circuit working  
c) Bushings for measuring relevant signals  
d) Piloting card for stepping motor with USB interface

**Technical Specification :**

Control Software  
USB cable  
Required services:  
Electrical supply: 220 V single-phase, 50/60 Hz.  
Dimensions:  
Power supply unit: 450 x 350 x 150 h mm  
Mechanical unit: 370 x 180 x 160 h mm  
Weight: 10 kg.



**Laboratory Instrument India**