



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stepping Motor Control Module

Product Code :
LBNY-0005-12600027



www.laboratoryinstrumentindia.com

Description :

Stepping motor, 200 step/rev, max torque = 60 N· cm
No. 2 limit micro switches
Graduated disc, 0 360
Encoder (200 pulse/rev)
Gear motor 90/1 rev.
Electric supply box consisting of:
a) Ammeter for measuring the motor absorbing current
b) Diagram on front panel showing the circuit working
c) Bushings for measuring relevant signals
d) Piloting card for stepping motor with USB interface

Technical Specification :

Control Software
USB cable
Required services:
Electrical supply: 220 V single-phase, 50/60 Hz.
Dimensions:
Power supply unit: 450 x 350 x 150 h mm
Mechanical unit: 370 x 180 x 160 h mm
Weight: 10 kg.



Laboratory Instrument India