



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hydraulic Pick And Place Modular Robot

**Product Code :**  
LBNY-0005-12600030



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It allows to perform educational exercises of programming, control and regulation on the axis movement. The operation phases correspond to the industrial techniques of loading and unloading used on the flexible lines. The Hydraulic Pick And Place Modular Robot robot has been designed with the proportional and on/off hydraulic technologies, electronic control. The unit is supplied together with a complete handbook.

**Technical Specification :**

**Technical Features:**

Max transportable load: 2.5 Kg  
Wrist rotation: 180°  
A axis rotation: 120°  
Working pressure: max 40 bar  
Z axis linear stroke: 150 mm.  
Electric requirements: 220/240 V - 50 Hz.



**Laboratory Instrument India**

