



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Axis Control Unit

**Product Code :**  
LBNY-0005-12600032



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The electrical equipment supplies the electrical feeding to all the components of the system, allows the manual command of the single axis or their automatic control through programs of complex cycles. The system allows to simultaneously or independently control two axes and to face the most important basic problems concerning industrial automation and specifically robotics.

**Technical Specification :**

**Required Services:**

Electric Power Supply 220/240 V AC - 50/60 Hz

Linear and Rotary Assembly with Support

Dimensions: 960x360x440h mm

Weight: 16 kg

Electrical Equipment

Dimensions: 540x470x360h mm

Weight: 28 kg



**Laboratory Instrument India**

