



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Basic Pneumatics Trainer

**Product Code :**  
LBNY-0005-12600033



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The trainer is supplied with an exhaustive manual including a technical description of all the components operating and maintenance instructions and many exercises.

The unit Basic Pneumatics Trainer is the economic version of the Modular Pneumological Tutor.

The unit is made of painted steel plate and can be utilised on any support bench.

The unit requires an external pneumatic supply. A compressor is supplied as a optional.

The unit has been designed to study and experiment the fundamental pneumatics.

**Technical Specification :**

**Required Services:**

Pneumatic supply max: 10 bar

Dimensions: 250x400x600

Net weights 20 kg

**Composition and Description:**

The trainer consists of a top-desk panel carrying the following components, industrial type.



**Laboratory Instrument India**

