



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Basic Oil Power Hydraulics Trainer

**Product Code :**  
LBNY-0005-12600035



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The Basic Oil Power Hydraulics Trainer is supplied with an exhaustivemanual, including a technical description of all the components, operating andmaintenance instructions and manyexercises. The unit Basic Oil Power Hydraulics Trainer is the economic version of the Modular oilpower hydraulic trainer. The panelrequires an external oil power supply unit, supplied as a optional. The unit ismade of a panel carrying the industrial type components needed to study basi coil power hydraulics. The Basic Oil Power Hydraulics Trainer consists of a top desk panelcarrying the following components:

**Technical Specification :**

**Required Services:**

Electric power supply: single-phase 220V AC 50/60Hz.

Hydraulic supply ca 9 l/min, a max. pressure 30bar.

Low viscosity 3-4Å°E hydraulic oil.

Dimensions: 1200 x 1000 x 800 h mm.

Net weight: 80 kg.



**Laboratory Instrument India**

