



Email : sales@laboratoryinstrumentindia.com

Product Name :
Angular Speed Module

Product Code :
LBNY-0005-12600036



www.laboratoryinstrumentindia.com

Description :

The automatic data acquisition and analysis system for Personal Computer allows to take the measurements automatically plot on screen or print out and file the test data. The Angular Speed Module module allows to study two standard measuring systems for the angular speed of shafts, such as the generator and the phonic wheel.

Technical Specification :

Required services:

Absorbed power: 0.1 kW max
Electric supply: 230 V, 50/60 Hz, single phase
Dimensions: 470 x 450 x 150 h mm
Net weight: 6 kg.



Laboratory Instrument India