



Email : sales@laboratoryinstrumentindia.com

Product Name :
Magnetic Damping Accessory To Add Friction On Cart

Product Code :
LBNY-0004-780013



Description :

Magnetic Damping Accessory To Add Friction On Cart

Technical Specification :

Mechanical bench accessory 11:

The Magnetic Damping Accessory can be used with either a Dynamics Cart, Smart Cart. The steel bracket of the Accessory is attracted to the magnets built into the cart (or the Magnetic Bumper on the Smart Cart), and the bracket easily slides up and down to adjust the amount of damping.

The magnetic drag is caused by the induced currents in the aluminum track, and is proportional to the speed of the cart. This accessory will not work with plastic tracks!

Magnetic damping accessory to add friction on cart, supplied with:

Magnets

Keeper

Adequate bracket.



Laboratory Instrument India

