



Email : sales@laboratoryinstrumentindia.com

Product Name :
Modular Pneumological Tutor

Product Code :
LBNY-0005-12600037



www.laboratoryinstrumentindia.com

Description :

The Modular Pneumological Tutor has been designed to study and experiment the fundamental pneumatics the pneumological techniques and the cascade pneumatic circuits.

Support Panel

Support Panel is made of painted steel plate and can be utilized on any support bench in sloping or horizontal position.

Technical Specification :

The Modular Pneumological Tutor is composed of the following parts:

Support panel

Module LPN for the study of the pneumatic logic

Module EPN for the study of the electro pneumatics

Module BPN for the study of the fundamental pneumatics.



Laboratory Instrument India

