



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Modular Oil Power Hydraulic Trainer

**Product Code :**  
LBNY-0005-12600038



#### **Description :**

It is suitable for technical and professional training at different levels. It consists of a trolley mounted structure with a hydraulic power supply unit and two panels, one for the basic oil power components and one for On/Off and proportional electrical control. The panels can be used independently, using an external hydraulic power supply with suitable characteristics. Modular Oil Power Hydraulic Trainer is a modular system specially designed to test out the techniques for the transfer of fluid power and electro hydraulic control techniques. The two panels are functionally independent, but can be coupled so as to form a double-faced structure, which is installed on the trolley mounted frame.

The trainer is supplied with an exhaustive manual, including a technical description of all the components, operating and maintenance instructions and many exercises. The Modular Oil Power Hydraulic Trainer can be easily interfaced with PLC or another programmable logic controller available to the user in order to program complex operating cycles.

#### **Technical Specification :**

Dimensions: 1200 x 1000 x 1600 h mm.

Net weight: 400 kg.

Gross weight: 610 kg.

---

# Laboratory Instrument India

