



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Transparent Oil Power Components Study Unit

**Product Code :**  
LBNY-0005-12600039



www.laboratoryinstrumentindia.com

**Description :**

The Transparent Oil Power Components Study Unit is composed of:

Overhead projector.

Set of transparent oil power components

Set of technical cards

Set of colour slides

Trolley-mounted furniture with support plane for the components storage, a drawer for tools, a support for projector

Over the furniture there is fitted an hydraulic power generator with the following characteristics:

Working pressure: 10 bar

Pressure limiting device

Manometer

Aluminium tank: 12 l

Single-phase electric engine 0.5 HP - 220 V single-phase 50 Hz

Air inlet valve in the circuit to evidence the flow.

**Technical Specification :**

**Power requirements:**

Electric power: 1kW

Main supply: 220 V - 50 Hz single-phase

Dimensions: 600 x 900 x 1000 h mm

Net weight: 60 kg.



**Laboratory Instrument India**