



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pneumatic Components Kit

Product Code :
LBNY-0005-12600040



www.laboratoryinstrumentindia.com

Description :

The sectioning procedure has been carefully studied to show the operating principles thus becoming a very important educational aid. Provides the components under sectioned version to visualize the technologies used their construction. You can find the same components fitted on our trainers for the study of pneumologic. Industrial type pneumatic components which are container in a suitable suitcase and are usually utilized in the industrial automation control circuits.

Technical Specification :

Weight and dimensions:

Suitcase dimensions: 450x320x120 mm

Net weight: 6 kg.



Laboratory Instrument India