



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electro Hydraulic Testing Station

Product Code :
LBNY-0005-12600043



Description :

Electro Hydraulic Testing Station learning system features a hydraulic press module, part clamp module, hydraulic power unit, and more. Among other demonstrations, the clamp module will pick up a test valve and the hydraulic press module will hold a test manifold against the valve body. Electro Hydraulic Testing Station allows learners to gain valuable product testing skills used in automated processes by studying topics like station operation and adjustment, module sequencing, and station sequencing. Learners will use these and other components to practice operating, programming, and adjusting real-world mechatronics equipment. Electro Hydraulic Testing Station learning system will allow learners to practice and study how products are tested on an automated line, how these skills are integrated within a larger automated process, and an example of how hydraulics are utilized on an automated line. Uses components that learners will find on-the-job in order to give the best opportunity to build confidence and industrial competencies.

Features:

Allows learners to practice and study how products are tested on an automated line, how these skills are integrated within a larger automated process, and an example of how hydraulics are utilized on an automated line

Allows learners to gain valuable product testing skills used in automated processes by studying topics like station operation and adjustment, module sequencing, and station sequencing

Includes a hydraulic press module with horizontal hydraulic cylinder and electro-hydraulic valve; hydraulic power unit; pressure transducer system; part clamp module; an operator station with lockout.

Technical Specification :

Electro Hydraulic Testing Station



Laboratory Instrument India