



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mechanical Oscillator to deliver signal for Harmonic Motion Experiments

**Product Code :**  
LBNY-0004-780014



**Description :**

Mechanical Oscillator to deliver signal for Harmonic Motion Experiments

**Technical Specification :**

Mechanical bench accessory 12:

The Mechanical Oscillator/Driver delivers repeatable, low-frequency, high-force sinusoidal motion for harmonic motion experiments. Shown with the Smart Cart Standard Dynamics System, it also works with other Dynamics Systems having a metal track.

Mechanical oscillator to deliver signal for harmonic motion experiments.

Amplitude with 12 cm

Mounted to Dynamic Track

Signal with at least of 10 VDC

Frequency of 3Hz

Current of 0.3 A



**Laboratory Instrument India**

