



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Power Supply

**Product Code :**  
LBNY-0005-12600054



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It has a voltage and current regulation and is therefore short circuit proof; it can supply DC 24V power from an additional pair of terminals. The Power Supply unit provides and monitors the power necessary for the system operation with safety extra low voltage: the integrated choke prevents the data telegrams from short circuiting on the bus line.

**Technical Specification :**

Power supply: 120 $\hat{A}$ -230 Vac, 5060 Hz  
Short circuit current: limited to 1.5 A  
Output voltage (Bus): 29 Vdc  
Output voltage: 29 Vdc



**Laboratory Instrument India**