



Email : sales@laboratoryinstrumentindia.com

Product Name :
Presence Detector And Brightness Sensor

Product Code :
LBNY-0005-12600055



Description :

The detector signal can be analyzed via two separate communication channels, termed motion detector and presence detector. Another main application is the automatic control of the lighting on an office workplace. Presence Detector And Brightness Sensor is a presence/motion detector with integrated constant light level control. It contains an independent light sensor with integrated 2nd level light control (switching) and constant light level control (dimming). Presence Detector And Brightness Sensor communicates via actuators or other devices.

Technical Specification :

Range: horizontal 360^o, vertical approx. 105^o, 288 sectors
Brightness measurement type: contrast
Power supply via bus line
Presence detection type: passive infrared (PIR)
Range: 201000 Lux
Current drain: approx. 10 mA
Different installation heights (from 2.5m to 5.0m).



Laboratory Instrument India

