



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Binary Output

**Product Code :**  
LBNY-0005-12600056



**Description :**

Each of the outputs can be assigned various tasks depending on the application program used: i.e. the Binary Output consists of the device (HW) and its application program (SW). The Binary Output is a device that can switch two separate groups of electric devices via its two outputs.

**Technical Specification :**

Power supply: Bus line  
Rated voltage: 230 Vac, 47-63 Hz  
Outputs: 2 volt free contacts  
Switching power: 1 kW (incandescent lamp) or 500 W (fluorescent lamp)  
Switching current at 24 Vdc: 10 A (resistive load), 4 A (inductive load).



**Laboratory Instrument India**