



Email : sales@laboratoryinstrumentindia.com

Product Name :
Binary Output

Product Code :
LBNY-0005-12600056



Description :

Each of the outputs can be assigned various tasks depending on the application program used: i.e. the Binary Output consists of the device (HW) and its application program (SW). The Binary Output is a device that can switch two separate groups of electric devices via its two outputs.

Technical Specification :

Power supply: Bus line
Rated voltage: 230 Vac, 47 \hat{A} -63 Hz
Outputs: 2 volt free contacts
Switching power: 1 kW (incandescent lamp) or 500 W (fluorescent lamp)
Switching current at 24 Vdc: 10 A (resistive load), 4 A (inductive load).



Laboratory Instrument India