



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Dynamics Track Mount

**Product Code :**  
LBNY-0004-780015



**Description :**

Dynamics Track Mount

**Technical Specification :**

Mechanical bench accessory 13:

Dynamics track mount:

Used to mount the rotary motion sensor and allows bi-directional movement on the track.

Supplied with a side rod to be fixed on the rotary motion sensor.

The Dynamics Track Mount is used to mount the Rotary Motion Sensor to the Dynamics Track, allowing it to act as a high resolution, bi-directional Smart Pulley.

When used in conjunction with the Dynamics Track Mount, the Track String Adapter allows a Rotary Motion Sensor to continuously monitor the Dynamics Cart position.

A loop of string wraps around the Rotary Motion Sensor pulley and the ball-bearing pulley, and then it attaches to the cart via a special clip (included).



**Laboratory Instrument India**

