



Email : sales@laboratoryinstrumentindia.com

Product Name :
Industrial Kneader Machine Training System

Product Code :
LBNY-0005-12600062



Description :

The kneading machine in this system is based on the machines that are commonly used in the food industry. Industrial kneading machines are widely used in food, medicine, rubber, paper and other industries. With reference to the actual teaching needs, the system includes the PLC control, variable frequency control, temperature control and touch screen control that are commonly used in industry, and covers the experimental trainings in sensor detection technology, food process control, power distribution technology and other subjects. They can realize stirring for a variety of raw materials and can also achieve heating at the same time.

General Features:

The trainer is composed of three independent modules:
The kneading machine,
The main electrical cabinet,
The control cabinet.

Technical Specification :

Industrial Kneader Machine Training System



Laboratory Instrument India

