



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Two Zone Parking Garage Simulator

**Product Code :**  
LBNY-0005-12600069



**Description :**

The automatic control is performed through a PLC, whose program manages the displays signaling vacant and full places for both zones and the opening and closing of the in and out barriers. Two Zone Parking Garage Simulator represents a two-zone parking garage where it is possible to realistically simulate the sequences that a car driver has to perform when he wants to use an automated parking garage. The number of cars in the parking garage is displayed through LED.

**Technical Specification :**

Two Zone Parking Garage Simulator



**Laboratory Instrument India**