



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Dead Time Element

**Product Code :**  
LBNY-0005-12600076



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Dead Time Element allows the insertion of an adjustable real dead time in those processes which are characterized by it.

**Technical features:**

Power supply: +15 V ; 0 V ; 15 V

Dead time  $T_t = 10 \text{ ms} \dots 100 \text{ ms} / 100 \text{ ms} \dots 1 \text{ s}$

Potentiometer fine setting.

Led indicators of over range.

Signal voltage range: •10V, ..., +10V

Proportional coefficient of the module  $KS = 1$

Three position switch for coarse setting and exclusion of the dead time.

**Technical Specification :**

Dead Time Element



**Laboratory Instrument India**

