



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dead Time Element

Product Code :
LBNY-0005-12600076



www.laboratoryinstrumentindia.com

Description :

Dead Time Element allows the insertion of an adjustable real dead time in those processes which are characterized by it.

Technical features:

Power supply: +15 V ; 0 V ; 15 V

Dead time $T_t = 10 \text{ ms} \dots 100 \text{ ms} / 100 \text{ ms} \dots 1 \text{ s}$

Potentiometer fine setting.

Led indicators of over range.

Signal voltage range: •10V, ..., +10V

Proportional coefficient of the module $KS = 1$

Three position switch for coarse setting and exclusion of the dead time.

Technical Specification :

Dead Time Element



Laboratory Instrument India

