



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mass and Hanger Set

**Product Code :**  
LBNY-0004-780017



**Description :**

Mass and Hanger Set

**Technical Specification :**

Mechanical bench accessory 15:

Mass and hanger set comprise:

4 mass hangers (5g each)

6 aluminium masses (3x10g + 3x5g)

12 brass masses (3x100g + 3x50g + 6x20g)

9 plastic masses (3x2g + 3x1g + 3x0.5g)

Supplied in a moulded plastic case



**Laboratory Instrument India**