



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Filling Tank

**Product Code :**  
LBNY-0005-12600079



**Description :**

Slide valves allow you to change the configuration of the controlled system. Filling Tank used to demonstrate the characteristics of the controller system. Also includes a set of 10 transparent sheets, with a graduated horizontal scale and with a vertical one, for a manual recording of the characteristics of the controlled system, suitable to the use of markers soluble in water. Includes plug in devices for the sensors used to measure the level of the liquid.

**Technical Specification :**

**Technical features:**

Tank capacity: 1



**Laboratory Instrument India**