



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sample and Hold Element

**Product Code :**  
LBNY-0005-12600085



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The sampling frequency can be provided by the generator which is integrated in the module or by an external signal. Used to discontinuously sample the behaviour of a continuous control on a process.

**Technical Specification :**

**Technical features:**

Sampling frequency: 0,2 .... 20 Hz

Signal voltage range: 10 V, ..., +10 V

Power supply: +15 V ; 0 V ; 15 V.



**Laboratory Instrument India**