




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Digital Input Simulator

**Product Code :**  
LBNY-0005-12600086



**Description :**

The board includes 32 switches on 4 columns for PLC input enabling, with a LED signaling the on/off state of the individual switch and the possibility of marking the references on special tags that can be cancelled; moreover, 4 outputs at 24 Vdc/1 A are available, internally powered. The simulator allows studying and analyzing the programming techniques through automatisms which are freely and imaginatively implemented. Moreover, it allows the PLC potentialities to be pointed out. Power supply: single-phase from mains. The connection to the PLC can be carried out through terminals or connector.

**Technical Specification :**

Digital Input Simulator



**Laboratory Instrument India**