



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Elastic Bumper

**Product Code :**  
LBNY-0004-780018



**Description :**

Elastic Bumper

**Technical Specification :**

Mechanical bench accessory 16:

The Elastic Bumper protects the Motion Sensor from the carts, but doesn't interfere with the ultrasonic pulse.

Elastic bumper:

Set of 2 brackets

10-meter at least of elastic material



**Laboratory Instrument India**