



Email : sales@laboratoryinstrumentindia.com

Product Name :
Derivative Action Element

Product Code :
LBNY-0005-12600087



www.laboratoryinstrumentindia.com

Description :

Derivative action controller suitable for closed loop continuous control systems.

Technical Specification :

Technical features:

Power supply: +15 V ; 0 V ; 15 V

Three position switch coarse setting.

Potentiometer fine setting.

Signal voltage range: 10V, ..., +10V

Led indicator of over range

Switch for inclusion/exclusion of the derivative action.

Coefficient of the derivative action $KD = 2 \text{ ms} \dots 2 \text{ s}$.



Laboratory Instrument India