



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Derivative Action Element

**Product Code :**  
LBNY-0005-12600087



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Derivative action controller suitable for closed loop continuous control systems.

**Technical Specification :**

**Technical features:**

Power supply: +15 V ; 0 V ; 15 V

Three position switch coarse setting.

Potentiometer fine setting.

Signal voltage range: 10V, ..., +10V

Led indicator of over range

Switch for inclusion/exclusion of the derivative action.

Coefficient of the derivative action  $KD = 2 \text{ ms} \dots 2 \text{ s}$ .



**Laboratory Instrument India**