



Email : sales@laboratoryinstrumentindia.com

Product Name :
Integral Action Element

Product Code :
LBNY-0005-12600088



Description :

Integral action controller suitable for closed loop continuous control systems.

Technical Specification :

Technical features:

- Three position switch coarse setting.
- Integral action reset input.
- Potentiometer fine setting.
- Power supply: +15 V ; 0 V ; 15 V
- Coefficient of the integral action $K_I = 0.1 \dots 100$
- Switch for inclusion/exclusion of the integral action.
- Led indicator of over range.
- Signal voltage range: 10V, ..., +10V.



Laboratory Instrument India