



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
P Controller

**Product Code :**  
LBNY-0005-12600089



**Description :**

Proportional action controller suitable for the closed loop continuous control systems.

**Technical Specification :**

**Technical features:**

Three position switch coarse setting.  
Led indicator of over range.  
Potentiometer fine setting.  
Power supply: +15 V ; 0 V ;15 V  
Proportional gain  $K_p = 0 \dots 100$   
Signal voltage range: 10V, .... , +10V.



**Laboratory Instrument India**