



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Gain And Offset Adjust

**Product Code :**  
LBNY-0005-12600091



**Description :**

Gain And Off set Adjust allows the proportional adaptation of the external signals to the normal voltages used in the automatic control systems.

**Technical Specification :**

**Technical features:**

Offset voltages that can be connected: 10 V .... +10 V

Potentiometer fine setting.

Coarse setting through rotary switches.

Power supply: +15 V ; 0 V ; 15 V

Voltage interval of the input signal:  $\hat{A}$ 50 V, ..., +50 V

Attenuation of the pulse signals.

Time constants: 0,1 .... 10 ms ; 10 .... 100 ms

Adjustable level through the setting of the gain: 0 ....1, 0 .... 10, 0 .... 100.



**Laboratory Instrument India**

