



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Matching Amplifier

**Product Code :**  
LBNY-0005-12600093



**Description :**

Amplifier as matching element between signal voltage levels and standard voltages used in automatic control systems.

**Technical Specification :**

**Technical features:**

Power supply: 15 V/0 V/+15 V

Coarse and fine gain setting: 01/010/0100

Input signal range U :50 V . . . + 50 V

Off set voltage: 10 V . . . +10 V

Connectable low pass filter with coarse and fine time constant setting: 0/ 110 ms / 10100 ms Connectable output.



**Laboratory Instrument India**