




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Programmable Logic Controller 26 In/22 Out

**Product Code :**  
LBNY-0005-1260095



**Description :**

Programmable Logic Controller 26 In/22 Out unit has technical and functional characteristics similar to the ones of the Programmable Logic Controller 26 In/22 Out, but it is composed of a CPU with 24 digital inputs, 2 analogue inputs and 22 relay outputs.

**Technical Specification :**

Programmable Logic Controller 26 In/22 Out



**Laboratory Instrument India**