




Email : sales@laboratoryinstrumentindia.com

Product Name :
Programmable Logic Controller 26 In/22 Out

Product Code :
LBNY-0005-1260095



Description :

Programmable Logic Controller 26 In/22 Out unit has technical and functional characteristics similar to the ones of the Programmable Logic Controller 26 In/22 Out, but it is composed of a CPU with 24 digital inputs, 2 analogue inputs and 22 relay outputs.

Technical Specification :

Programmable Logic Controller 26 In/22 Out



Laboratory Instrument India