



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Programmable Logic Controller 12 In/8 Out

**Product Code :**  
LBNY-0005-12600096



**Description :**

Simulation of the inputs through switches and externally through terminals. The programming software is also supplied to allow programming the PLC from the computer, if so preferred. Easy to be programmed from its own panel, without computer. It is composed (including the extension module) of 12 digital inputs and 8 relay outputs.

**Technical Specification :**

Programmable Logic Controller 12 In/8 Out



**Laboratory Instrument India**