



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Multimedia Learning Software Piston Compressor

**Product Code :**  
LBNY-0005-12700011



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Multimedia learning software as an accessory to assembly projects  
Simple operation with user interfaces similar to MS Windows  
Editable drawings and parts lists  
Access to parts lists and drawings.  
Password-protected configuration files.

**Technical Specification :**

**Learning Objectives / Experiments:**

Understanding engineering drawings, exercises with AutoCAD files  
Virtual assembly and disassembly  
Explanation of functions  
Printout of parts lists  
Computer-aided calculations  
Dealing with design aspects  
Planning the assembly sequence  
Familiarisation with single parts and modules  
Use of Technical English on a real project.

---

# Laboratory Instrument India

