



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multimedia Learning Software Piston Compressor

Product Code :
LBNY-0005-12700011



www.laboratoryinstrumentindia.com

Description :

Multimedia learning software as an accessory to assembly projects
Simple operation with user interfaces similar to MS Windows
Editable drawings and parts lists
Access to parts lists and drawings.
Password-protected configuration files.

Technical Specification :

Learning Objectives / Experiments:

Understanding engineering drawings, exercises with AutoCAD files
Virtual assembly and disassembly
Explanation of functions
Printout of parts lists
Computer-aided calculations
Dealing with design aspects
Planning the assembly sequence
Familiarisation with single parts and modules
Use of Technical English on a real project.

Laboratory Instrument India

