



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laminar Flow Reactor

Product Code :
LBNY-0005-12700012



Description :

Laminar flow reactor for connection to supply unit
Transparent double jacket from pmma for cooling and heating with
Sensor for measuring the conductivity and temperature via
Air vessel for damping of pulsation
Straight glass tube with laminar flow
Temperature control in the reactor via
T-piece with nozzle for mixing the reactants
Special inlet for reducing the inlet length.

Technical Specification :

Laminar flow reactor:
Material: glass
Height: 600mm
Inside diameter: 15mm

Water Bath
Material: PMMA
Capacity: approx. 0,45l
Inside diameter: 45mm.



Laboratory Instrument India