



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Discontinuous Stirred Tank Reactor

**Product Code :**  
LBNY-0005-12700013



**Description :**

Chambered bottom made of stainless steel as heat exchanger for connection to Sensor for measuring the conductivity and temperature via  
Reactor with stirrer  
Temperature control in the reactor via  
Discontinuous stirred tank reactor for connection to supply unit.

**Technical Specification :**

Reactor  
Capacity: 750ml  
Speed stirrer: 330min<sup>-1</sup>  
Outer diameter: 110mm  
Inside diameter: 100mm  
Height: 140mm.



**Laboratory Instrument India**

