



Email : sales@laboratoryinstrumentindia.com

Product Name :
Discontinuous Stirred Tank Reactor

Product Code :
LBNY-0005-12700013



Description :

Chambered bottom made of stainless steel as heat exchanger for connection to Sensor for measuring the conductivity and temperature via
Reactor with stirrer
Temperature control in the reactor via
Discontinuous stirred tank reactor for connection to supply unit.

Technical Specification :

Reactor
Capacity: 750ml
Speed stirrer: 330min⁻¹
Outer diameter: 110mm
Inside diameter: 100mm
Height: 140mm.



Laboratory Instrument India

